

How to install and start the US-Smart device on your tablet and Smartphone.

1 - You should pair the US-Smart on your tablet:

- I. To do this, go to the Bluetooth parameters on your device.
- II. Switch Bluetooth on and switch the US-Smart on, with the Power Switch:
(The Led of US-Smart becomes green)
- III. Start a device search and select "IOIO (BF:EB)" device. It is the US-Smart.
A password will be requested, enter: 4545
- IV. And that's it, the US-Smart is paired.

2 - For install the application:

- I. Go to our website: <http://www.lecoeur-electronique.com/>
- II. Click on the US-Smart picture or on the tab "Ultrasound Products" and "New device" for arrived on the US-Smart page
- III. On the right of the page, you will find all downloadable documents and applications.
- IV. Download it, install it and that's it.

If you need help or have any questions, don't hesitate to contact us at:

info@lecoeur-electronique.com