




Android Ultrasound Device

US-ANDROID

Datasheet V0 01 2020







US-ANDROID is our latest ultrasound device, it transforms your phone or your tablet in a high performances UT instrument.

US-ANDROID is our new generation ultrasound devices with a single channel to transmit and receive ultrasonic waves. Its very small size and its advanced technology allow to have a unique product for more applications like medical ultrasound imaging, the NDT and also for the research and university. The transmitter can generate pulses with a voltage level and a width programmed by the user. A low noise preamplifier combined to a VGA gives a gain range between 0 and 80 dB, a DAC curve is also available. A 12 bits analog digital converter with a sampling frequency of 80 MHz is used to digitize ultrasound signals. The device has 2 working modes: Transmission or Reflection. The power supply is delivered by the USB connector of the android device.





TECHNICAL SPECIFICATIONS

- Android OS (V4.00 ->)
- 4/6 Hours Autonomy (tested with Samsung S10 device)
- USB connection
- Tactile GUI
- Ultralow noise preamplifier: 0.74 nV/VHz
- 6 dB bandwidth: 540 KHz to 18 MHz
- High voltage transmitting pulses (10 to 230 Volts)
- 50 Ω load drive
- Single or double crystal operation (software selectable)
- Digitizer 12 bits at 80 MSPS
- Programmable gain: 0 to 80 dB
- DAC curve
- TGC curve
- USB powered
- Very small size: 70 x 36 mm







PICTURES

Typical application (android device/ US-ANDROID / probe)










les ultrasons






Main Menu








TGC









Other SC...







LECOEUR ELECTRONIQUE				Wednesda	iy, January 8	,2020 01:49
LE GAINTE 1						
K≓ SCALE ^{0.8 %}						
lī <u>-</u> GATES ₀₁						
TRANSMITTER						
ROBE						
MT TGC CURVE						
k⊉ CAD ^{std}						
MEASUREMENT UNIT / VELC	CITY					
₩ SIGNAL ANALYSIS TOOLS						
SAVE/RECALL SETUP						
UPTIONS						
P QUITNDROID < 0 µS			8.95 µS/di	V		89.46 µS ->
	III	0	<			



les ultrasons



les ultrasons



